

MWi PHARM++
The best choice for TDM.

MWi PHARM++
The best choice for TDM.

Among all available Therapeutic Drug Monitoring software MWiPHARM++ has the largest database:

- 175 drugs with their pharmacokinetic properties
- 285 population models

MWi PHARM++
The best choice for TDM.

MwPharm++ is an efficient application for Therapeutic Drug Monitoring (TDM) to optimize pharmacotherapy by establishing a proper dosing regimen based on population PK parameters from an extensive drug database and individual physiological patient parameters

MWi PHARM++
The best choice for TDM.

20+ years
MWiPHARM is trusted in clinical practice for 20+ years

1991
Dutch Standard since 1991 by the Dutch Association of Hospital Pharmacists (NVZA)

MWi PHARM++
The best choice for TDM.

Administrator Medication

Observe Blood Samples

Analyze in MwiPharm++

Optimize Dosing Regime

MWi PHARM++
The best choice for TDM.

MWPHARM has been awarded in 2013 to be the world's best solution for TDM software

Ruth A. Gude & Thomas Sauer J. Wilmer. Monitoring therapeutic drug monitoring software: a review of available computer tools. Clinical Therapeutics. 2013 Jan;35(1):10-20. doi: 10.1007/s10280-012-0800-y.

MEDIWARE

MW PHARM++

The best choice for TDM.

Among all available Therapeutic Drug Monitoring software MWPHARM++ has the largest database:

- 175 drugs with their pharmacokinetic properties
- 285 population models

MwPharm++ is an efficient application for Therapeutic Drug Monitoring (TDM) to optimize pharmacotherapy by establishing a proper dosing regimen based on population PK parameters from an extensive drug database and individual physiological patient parameters

Obtain blood samples → Analyze in MwPharm++ → Administer/Verification → Optimize Dosing Regime

MW PHARM++
The best choice for TDM.

20+ MWPHARM is trusted in clinical practice for 20+ years

1991 Dutch Standard since 1991 by the Dutch Association of Hospital Pharmacists (NVZA)

NVZA

MW PHARM++
The best choice for TDM.


MWPHARM has been awarded in 2013 to be the world's best solution for TDM software

Rich A, Coglio C, Thome V, Bax J, Wilmar S. Optimizing therapeutic drug monitoring software: a review of available computer tools. *Clinical Pharmacology and Therapeutics*. 2014;95(1):10-16. doi:10.1002/cpt.100

MEDIWARE



MwPharm++ is an efficient application for Therapeutic Drug Monitoring (TDM) to optimize pharmacotherapy by establishing a proper dosing regimen based on population PK parameters from an extensive drug database and individual physiological patient parameters



System Requirements

Minimum Requirements

Processor: 1.2 GHz or higher
(32-bit (x86) or 64-bit (x64) processor)

RAM: 512 MB or more
(available for the application)

Harddisk space: 1 GB

Supported Operating Systems:
Windows Vista SP2, Windows 7,
Windows 8, Windows 8.1

Virtual desktop environments such as Citrix
are also possible after configuration





The best choice for TDM.

Among all available Therapeutic Drug Monitoring software
MWPHARM++ has the largest database:

175 drugs with their pharmacokinetic properties

285 population models.




The best choice for TDM.

20+
years

MWPHARM is trusted in clinical practice
for 20+ years

1991

Dutch Standard since 1991
by the Dutch Association of Hospital
Pharmacists (NVZA)





The best choice for TDM.



**MWPHARM has been awarded in 2013
to be the world's best solution for TDM software**

Fuchs A, Csajka C, Thoma Y, Buclin T, Widmer N.:Benchmarking therapeutic drug monitoring software: a review of available computer tools. Clinical Pharmacokinet. 2013 Jan; 52(1):9-22. doi: 10.1007/s40262-012-0020-y



The best choice for TDM.



Dutch Standard since
1991

NVZA - the Netherlands



CE-marking
Software conforms to
legal requirements

Europe



ISO 13485 Certified
quality management system for
medical devices

Europe



ISO 9001 Certified
quality management system for
organisations

Europe



The best choice for TDM.



MwPharm DOS has been casting its light for more than 20 years

MWPHARM++ has been greatly improved to have it rock solid for the next 20+ years. For this, a new foundation has been developed for MWPHARM++ called EDSIM++



The approach to TDM you are familiar with



The new foundation for MWPHARM++ with new PK/PD simulation capabilities




Log in

Additional users can be created under the User tab.

The screenshot shows a software window titled 'MwPharm++'. At the top, there is a menu bar with File, Model, Chart, Macro, Object, and Help. Below the menu is a toolbar with various icons. A search bar contains 'Name:' and 'Type:'. The main area is titled 'Model' and shows a list of users:

Name	Username	Password	Level
Administrator			
Super User			
Normal User			

At the bottom right of the window is an 'Insert' button.





The best choice for TDM.

The approach to TDM you are familiar with

"Ease of transition for existing users is one of the focus areas for MwPharm++"

While MwPharm++ is built on the new foundation Edsim++,

this version has the look & feel you are familiar with

(see for example, the Patient tab on the right),

to ease the transition from MwPharm DOS to MwPharm++.

Of course, shortcut keys (such as F12) are still present in MwPharm++.


The screenshot shows the MwPharm++ software interface with the 'Patient' tab selected. On the left, there is a form with fields for Patient Number (1234567890), Name and Initials (Johnson, D.), Date of Birth (18-10-1955), Sex (Female), Address, Postcode / Zipcode, City, Family Doctor, Requesting Physician, Ward, Room Number, Description, Medication Date (21-3-1990), Age (34 years), Last Medication (gentamicin), and Date of Change (11-10-2015). To the right of the form is a large list of patient records. The first few columns of the list are 'Dob' (Date of Birth), 'Number', and 'Name'. The list contains numerous entries, such as 11-11-1915, !C007, !CASUS 007, 13-9-1990, !C008, !CASUS 008, 11-11-1924, !C009, !CASUS 009, 11-11-1977, !C010, !CASUS 010, 26-4-1926, !N001, !NPEM 001, 27-12-1929, !N002, !NPEM 002, 26-4-1923, !N003, !NPEM 003, 22-11-1984, !N004, !NPEM 004, 22-3-1921, !N005, !NPEM 005, 31-10-1942, !N006, !NPEM 006, 15-4-1910, !N007, !NPEM 007, 30-8-1926, !N008, !NPEM 008, 11-5-1930, !N009, !NPEM 009, 27-12-1949, !N010, !NPEM 010, 19-7-1939, 0 EXAMPLE 1, Adult, 19-7-1939, 0 EXAMPLE 2, Dialysis, 19-7-1939, 0 EXAMPLE 3, CAPD, 14-7-1994, 0 EXAMPLE 4, Neonate, 19-7-1939, 0 HISTORY, Tour De Pharm, 1-1-1951, 102, TEST-D, 1-1-1935, 1234, TEST-E, 18-10-1955, 1234567890, Johnson, D., 1-1-1979, 1337, Henk Test, 1-1-1979, 13371337, Ingrid Test, 1-1-1958, 31, TEST-B, 1-1-1958, 51, TEST-A, 12-9-1948, 65, TEST-C, 1-1-1961, 8, DEMO REGIMENT, 1-1-1935, 88, DEMO DIALYSIS. At the bottom of the list are buttons for Standard, Load, Active, Delete, Update, and Insert. Below the list is a navigation bar with links like Login, Patient, Status, Case, History, Simulation, Fitting, Dosing, Users, Settings, Help, and About. The word 'Sim' is highlighted in the status bar at the bottom.

Status tab

After selecting a patient, you can fill in the physiological patient parameters in the Status tab.

Select for example: the Renal Function of the patient

You can also see a small summary of the Population model here



Manage cases of patients with ease

After selecting a patient:


Select one of the 175 drugs and one of the 285 population models

Click on *Insert* to create a Case

You can add as many cases to the same patient as you prefer.

In MwPharm DOS it was common to create an additional record of the same patient (and include date and/or patient number in the name of the record) for creating a new dosing regime.

The benefit of creating Cases is that the list of patients under Patient tab will contain only unique patients.







The best choice for TDM.

The approach to TDM you are familiar with

You can still deselect a row using the Space bar on the keyboard.

Hemodialysis can also be added.


You can find the familiar shortcut keys under the red triangles.



Simulation tab


Generate a simulation based on the Dosage regimen created in the History tab.

Create an advise with the Text tool




Simulation tab

Toggle check boxes to show other parameters such as Area Under Curve and Area Under Curve Cumulative.



Generate a Report (via R button).
This can create a PDF that can be saved, printed, e-mailed et cetera.




Fitting tab

Generate a Fitting

See the difference between
the unfitted simulation and
the fitted simulation.

The unfitted simulation



Fitting tab

3. Lastly, press Fit to generate a Fitting with Bayesian Parameters included.

1. Generate Bayesian Parameters

2. After generating the Bayesian Parameters, set type to Bayesian.

Standard Deviation (SD)
This is what the Bayesian Parameters will generate based on the population model.




Dosing tab

Generated optimized dosing regime (black) based on Fitting.

Dosing regime (red) based personal settings.
(in this case: 24 hour, 8 dosages)

Change dosing regime to preferred settings. In this case:
every 24 hours and 8 dosages.

Add the dosage regime with preferred settings to History
(will be added to a new row).



MW PHARM++

The best choice for TDM.

Among all available Therapeutic Drug Monitoring software MWPHARM++ has the largest database:

- 175 drugs with their pharmacokinetic properties
- 285 population models

MwPharm++ is an efficient application for Therapeutic Drug Monitoring (TDM) to optimize pharmacotherapy by establishing a proper dosing regimen based on population PK parameters from an extensive drug database and individual physiological patient parameters

Obtain blood samples → Analyze in MwPharm++ → Administer/Verification → Optimize Dosing Regime

MW PHARM++
The best choice for TDM.

20+ MWPHARM is trusted in clinical practice for 20+ years

1991 Dutch Standard since 1991 by the Dutch Association of Hospital Pharmacists (NVZA)

NVZA

MW PHARM++
The best choice for TDM.

MWPHARM has been awarded in 2013 to be the world's best solution for TDM software

Ricke A, Coglio C, Thome V, Bax J, Mihnev N. Optimizing therapeutic drug monitoring software: a review of available computer tools. *Clinical Pharmacology and Therapeutics*. 2013 Jan;93(1):14-21. doi:10.1002/cpt.100

MEDIWARE